


**PROPER SPINE ALIGNMENT
IS THE KEY TO THE PRESENTATION
AND SELECTING THE RIGHT SLEEP SYSTEM.**



The correct mattress depth will allow the shoulder and hip to sink deep enough into the mattress to stabilize the spine while minimizing pressure point pain on the shoulder and hip.

The shoulders and hips will be perpendicular and the spine straight or in a neutral position if the mattress is correct.

When in doubt - Checking Posture is 1st priority

Knowing their room temperature is 2nd



The **Martenique Supreme** has its 10-inch profile designed for smaller framed persons looking for a firmer support system that enhances comfort and responsiveness. It has 3 inches of 4 lb memory foam on a 7-inch high-density foam support core.

Supreme sleep surface does not harden in cooler room temperatures insuring better spinal support, minimizing pressure point pain and unrestricted blood circulation.

It is priced as an entry-level sleep system. The removable washable cover is allergen-resistant.

Competes against Tempur-pedic firmest 8-inch profile OriginalBed and AdvantageBed.



Martenique Deluxe 550 with its 10-inch profile combines the soft, responsive layer of our latest breakthrough advancement – ventilating a 5-inch layer of our 2nd generation 5.25 lb memory foam with a 5-inch high density foam base layer creating a soft contouring support.

The ventilated memory foam increases air circulation so you sleep cooler, more comfortable compared to other national memory foam brands.

It is designed to provide a softer feel plus the added responsiveness for easier turning, proper spinal alignment, minimize pressure point pain, no motion disturbance and unrestricted blood circulation.

The removable washable cover is allergen-resistant.

This mattress can accommodate people over 250+ lbs The mattress has proven to be a great choice for expecting mothers, fibromyalgia suffers, degenerative disc, bludging or herniated disc and rheumatoid arthritis.

Competes against Tempur-pedic newest CloudBed line and the DeluxeBed.



Martenique Solution Bed is our support mattress and is a 10-inch profile consisting of 4 inches of 4.5lb 2nd generation memory foam and a 6-inch 2.0lb polyurethane support core. The design of this mattress creates a faster recovery memory foam making it easier to turn on and for transferring limited mobility patients.

It is designed to provide added responsiveness for proper spinal alignment, minimize pressure point pain, no motion disturbance and unrestricted blood circulation.

The Outlast cover, designed by NASA, is designed to disperse body heat and is allergen-resistant and anti-dust mite.

The foam density combination creates our best support mattress for a number of spinal issues including sciatica, scoliosis, hip replacement, degenerative disc, bludging or herniated disc and rheumatoid arthritis.

Competes against Tempur-pedic newest CloudBed line, the Classic Bed, the DeluxeBed and the RhapsodyBed.



Martenique Whisper Escape is an 11-inch profile with 1 inch of 4 lb 2nd generation memory foam quilted into the cover with 4 inches of 4.5 lb memory foam supported by a 6-inch 2.0lb polyurethane support core. The design of the surface creates a faster recovery with less pressure point making it easy to turn on.

The mattress was designed to minimize pressure points while maintaining good spinal support.

The 2nd generation memory foam is not affected by room temperatures, which makes it a good solution for fibromyalgia sufferers who sleep in a colder bedroom.

The Outlast cover, designed by NASA, is designed to disperse body heat and is allergen-resistant and anti-dust mite.



GelLogic Majestic has a 14-inch profile designed to relieve chronic back pain. The fusion of form and function is created when the proprietary gel is integrated with the foam through a patented process. The unique pattern provides true orthopedic spinal support and pressure reduction at the same time.

Unlike visco-elastic foam, or more commonly known as memory foam, the gel pattern combined with our 2nd generation memory foam immediately responds to shifts in weight. This technology creates an unsurpassed comfort layer that allows for easier repositioning and turning while at the same time does not telegraph your partner's motion disturbing your sleep. The new technology creates a cooler sleep environment with its heat sink properties absorbs body heat instead of reflecting body heat insuring a better night's sleep.

Our uniquely formulated memory gel sleep surface conforms to every curve and angle of your body with exacting precision, providing more support, greater pressure relief and long lasting durability.

The gel/foam combination creates a fuller feeling support in the lumbar and spinal support for a number of spinal issues including sciatica, scoliosis, hip replacement, degenerative disc, bludging or herniated disc and rheumatoid arthritis.

The gel, foams and Outlast cover used to make this ergonomic sleep system is hypoallergenic, anti-bacterial and anti-dust mite.

iC 900



Comfortaire iC900 has a 15-inch profile and is designed to provide a more responsive adjustable air mattress compared to Sleep Number.

The iC900 components as shown below has the new edge-to-edge sleep surface.

The anti-fungus and anti-mold 3.5-inch air bladder is more responsive for better support compared to Sleep Number.

The overall design creates even more support and unsurpassed durability compared to Sleep Number

It is an excellent alternative to Sleep Number or any inner spring mattress.

▶ **Mattress Cover**

Natural bamboo cover, a highly absorbent mesh-like structure with microscopic gaps in the fiber to ensure the fabric always feels cool and fresh, with suede sides.

▶ **Topper Foam**

Plant based synthetic latex 2¾-inch convoluted topper foam.

▶ **Air Chambers**

3½-inch urethane alloy air cores with patented design.

▶ **Foam Perimeter Support**

1.6 density foam with fabric covers and Velcro fasteners. 6" high

▶ **Base Foam**

Plant based synthetic latex cushion with fabric covers.

▶ **Air Control System**

Whisperflo® 4200 air control system with 100 digitally controlled settings



Alzheimer – Supreme for people under 150lbs Deluxe 550 or Solution for over 150lbs. Possibly the Whisper Escape.

BACK

Chronic Back Pain – SolutionBed or GelLogic Majestic

Bludging or herniated disc – Deluxe 550, SolutionBed or GelLogic Majestic. Possibly the iC900. Supreme would be too hard and Whisper Escape too soft.

Degenerative disc- Deluxe 550, SolutionBed or GelLogic Majestic. Possibly the iC900. Supreme would be too hard and Whisper Escape too soft.

Cancer – Chemotherapy – Deluxe 550, SolutionBed or GelLogic
Possibly iC900.

Toppers – 3” memory foam or if the mattress is breaking down the GelLogic topper. Keep in mind for a topper to work it needs to be placed on a firm mattress. If the mattress is too soft it will not provide the needed support combination for the topper to work.

Fibromyalgia – Deluxe 550 or Whisper Escape
Deluxe 550 or Solution if there are spinal issues combined with Fibromyalgia

Hip replacement –Martenique Solution is the first choice. A great alternative is the GelLogic Majestic particularly for bigger people (over 225lbs). The Supreme could work for smaller frame patients but the Whisper Escape is too soft.

Limited Mobility –GelLogic Majestic or Solution.

Pregnancy - Deluxe 550, SolutionBed or GelLogic Majestic. Both beds provide the right support with comfort to help minimize symptoms. Possibly the iC900.

Rheumatoid arthritis - Deluxe 550, SolutionBed or GelLogic Majestic.
Both beds provide the right support with comfort to help minimize symptoms.

Possibly the iC900. Supreme would be too hard and Whisper Escape too soft.

Sciatica – All the beds help depending on the size of the person – the person's posture needs to be checked while lying down to determine which will work the best.

Small framed patients can consider using the Supreme. Larger patients can choose between Deluxe 550 or Martenique Solution or GelLogic Majestic.

Scoliosis – Deluxe 550 or Martenique Solution

Remember when in doubt - Checking Posture is 1st priority

Knowing their room temperature is 2nd